

Lecture presentation

Open Access

Single-molecule tracking of raft-based signal transduction: a system of digital signal transduction?

Akihiro Kusumi

Address: Membrane Mechanisms Project, ICORP-JST, Institute for Integrated Cell-Material Sciences (iCeMS); Research Center for Nano Regenerative Medical Engineering, Institute for Frontier Medical Sciences, Kyoto University, Kyoto, Japan

Email: Akihiro Kusumi - akusumi@frontier.kyoto-u.ac.jp

from Seventeenth Annual Computational Neuroscience Meeting: CNS*2008
Portland, OR, USA. 19–24 July 2008

Published: 11 July 2008

BMC Neuroscience 2008, 9(Suppl 1):L6 doi:10.1186/1471-2202-9-S1-L6

This abstract is available from: <http://www.biomedcentral.com/1471-2202/9/S1/L6>

© 2008 Kusumi; licensee BioMed Central Ltd.

The plasma membranes of neuronal cells contain high concentrations of glycosphingolipids and glycosylphosphatidylinositol-anchored receptors (GPI-ARs), as well as cholesterol, which suggests important roles played by hypothetical microdomains, called raft domains in these membranes. Using simultaneous two-color single-molecule tracking of GPI-ARs, as well as intracellular lipid-anchored signaling molecules, G α i, Lyn, and PLC γ , we have obtained results showing that the plasma membrane is poised for assembly of these molecules, upon the external stimulation that initiates oligomerization of 3–9 GPI-AR molecules.

The receptor-cluster-induced, cholesterol-dependent assembly, termed receptor-cluster raft (RCR), works as a platform for the signal transduction of GPI-AR. G α i2 and Lyn (GFP conjugates) are recruited to RCRs frequently, but transiently (100–200 ms), based on protein-protein and lipid-lipid (raft) interactions. G α i2 binding to and its subsequent activation of Lyn are likely to take place within the same RCR, resulting in actin-dependent temporary immobilization (0.57-s lifetime, called Stimulation-induced Temporary Arrest of Lateral diffusion or STALL, every 1.3 s), inducing the temporary (250 ms) recruitment of PLC γ 2, for IP $_3$ production. Therefore, the RCR in STALL is a key, albeit transient, platform for transducing the extracellular GPI-AR signal to the intracellular IP $_3$ -Ca $^{2+}$ signal, via PLC γ 2 recruitment.

The bulk activation of IP $_3$ -Ca $^{2+}$ signaling and Src-family kinases persists over several minutes to several 10 s of minutes. Meanwhile, single-molecule events, such as

STALL and the recruitment of PLC γ 2, G α i2, and Lyn to RCR, lasted only for a fraction of a second. Namely, individual single-molecule events may occur like a digital pulse, and the bulk analogue-type activation of signaling molecules may be the result of superposition of these pulse-like signals. In this sense, the basic signaling mechanism in the raft-based signaling system could be called digital or frequency-modulated.